

Individual cell-based virological assays require a maximum of experience and know-how. After different viral infections in vitro, the cytopathic effect (CPE) occurs, which changes the shape and appearance of the infected cells. In the early infection, only few cells show a CPE, which makes it difficult i.e. to estimate the success of an infection experiment. In later stages, quantification and qualification of CPE showing cells can be helpful to assess the right timepoint for different experimental approaches.

With the help of the AID multiSpot (Fig. 4), a combination of an automated inverted microscope and a fluorescence EliSpot reader, this CPE can be monitored and quantified without cellular staining procedures and without interrupting any experimental approach directly in plates with liquid like medium.




Fig.1 Examples of DMSO treated (B,D) or untreated cells (A,C), comparison of transmitted light (A,B) and phase contrast (C,D) in 10x magnification




Fig.2 Dynamic measurement of vaccinia infected cells on (A) day 9 (CPE 0), (B) day 10 (CPE 1) and (C) day 11 (CPE 4)




Fig.3 (A) viral plaque assay, well view, with counting dialogue example (red lines). (B) infected cells stained with DAPI and (C) GFP

Summarizing, this approach for the AID multiSpot enables to quantify and visualize the CPE in virus infected or compound treated cells and to get a complete overview of the whole wells that no important event is missed. This can be done in order to relate different treated cells showing a CPE and to make results between different experiments, concentrations or timepoints comparable.



Fig.4 The AID multiSpot with different examples

